

```

# =====
# NOTEBOOK 2: MODEL TRAINING (v17 - FIXED: Acc + Early Stopping + Plots)
# =====

import os
import re
import json
import time
import random
import zipfile
import numpy as np
import pandas as pd
import cv2
import torch
import torch.nn as nn
import torchvision.transforms as T
import timm
from pathlib import Path
from PIL import Image
from dataclasses import dataclass
from torch.utils.data import Dataset, DataLoader, WeightedRandomSampler
from sklearn.metrics import f1_score, classification_report, confusion_matrix, roc_auc_score
from sklearn.model_selection import GroupShuffleSplit
import matplotlib.pyplot as plt
import seaborn as sns
from google.colab import drive

# 1. Setup
drive.mount('/content/drive')

# --- CONFIGURATION ---
@dataclass
class CFG:
    ZIP_PATH = "/content/drive/MyDrive/MLProject/Cleaned_Dataset_Processed.zip"
    EXTRACT_T0 = "/content/via_dataset_ready"

    TRAIN_DIR = os.path.join(EXTRACT_T0, "train")
    TEST_DIR = os.path.join(EXTRACT_T0, "test")
    CSV_PATH = "/content/drive/MyDrive/MLProject/labels_combine.csv"

    NUM_CLASSES: int = 3
    IMAGE_SIZE: int = 384
    BATCH_SIZE: int = 32
    EPOCHS: int = 25 # Increased limit, Early Stopping handles the rest

    LR: float = 1e-5

    SEED: int = 42
    MODEL_NAME: str = "efficientnet_b0"
    OUT_DIR: str = "runs"

    # Regularization
    DROPOUT_RATE: float = 0.5
    DROP_PATH_RATE: float = 0.2
    PATIENCE: int = 5 # Stop if no improvement for 5 epochs

    ZOOM_FACTOR: float = 0.65

cfg = CFG()

def set_seed(seed=42):
    random.seed(seed); np.random.seed(seed)
    torch.manual_seed(seed); torch.cuda.manual_seed_all(seed)
    torch.backends.cudnn.deterministic = True

set_seed(cfg.SEED)
device = torch.device("cuda" if torch.cuda.is_available() else "cpu")
print(f"Device: {device}")

# 2. Extract Data
if os.path.exists(cfg.EXTRACT_T0):
    import shutil
    shutil.rmtree(cfg.EXTRACT_T0)

with zipfile.ZipFile(cfg.ZIP_PATH, 'r') as z:
    z.extractall(cfg.EXTRACT_T0)

# 3. Data Preparation
def prepare_dataframes():
    master_df = pd.read_csv(cfg.CSV_PATH)
    label_map = dict(zip(master_df['filename'], master_df['label']))

    def scan_folder(folder_path):
        data = []
        if not os.path.exists(folder_path): return pd.DataFrame(data)
        for f in os.listdir(folder_path):
            if not f.lower().endswith('.jpg'): continue

            clean_name = re.sub(r'_aug.*', '', Path(f).stem)
            original_key = f"{clean_name}.jpg"

            # Robust Patient ID extraction
            if "sample" in clean_name:
                patient_id = clean_name.split('_')[0]
            else:
                match = re.search(r"([0-9]+[A-Z]+[A-Z]+[0-9]+)", clean_name)
                patient_id = match.group(0) if match else clean_name

    def scan_folder(folder_path):
        data = []
        if not os.path.exists(folder_path): return pd.DataFrame(data)
        for f in os.listdir(folder_path):
            if not f.lower().endswith('.jpg'): continue

            clean_name = re.sub(r'_aug.*', '', Path(f).stem)
            original_key = f"{clean_name}.jpg"

            # Robust Patient ID extraction
            if "sample" in clean_name:
                patient_id = clean_name.split('_')[0]
            else:
                match = re.search(r"([0-9]+[A-Z]+[A-Z]+[0-9]+)", clean_name)
                patient_id = match.group(0) if match else clean_name

```

```

        if original_key in label_map:
            data.append({
                "path": os.path.join(folder_path, f),
                "label": label_map[original_key],
                "patient_id": patient_id
            })
    return pd.DataFrame(data)

train_full_df = scan_folder(cfg.TRAIN_DIR)
test_df = scan_folder(cfg.TEST_DIR)

mapping = {"Normal": 0, "Precancerous": 1, "Cancerous": 2}

train_full_df['class_num'] = train_full_df['label'].apply(lambda x: mapping.get(str(x).capitalize(), mapping.get(str(x), -1)))
test_df['class_num'] = test_df['label'].apply(lambda x: mapping.get(str(x).capitalize(), mapping.get(str(x), -1)))

train_full_df = train_full_df[train_full_df['class_num'] != -1]
test_df = test_df[test_df['class_num'] != -1]

# Validation Split
gss = GroupShuffleSplit(n_splits=1, test_size=0.15, random_state=cfg.SEED)
train_idx, val_idx = next(gss.split(train_full_df, groups=train_full_df['patient_id']))

return train_full_df.iloc[train_idx].reset_index(drop=True), train_full_df.iloc[val_idx].reset_index(drop=True), test_df

# 4. Transforms (CLAHE + Zoom)
class ApplyCLAHE:
    def __init__(self, clip_limit=2.0): self.clahe = cv2.createCLAHE(clipLimit=clip_limit, tileGridSize=(8, 8))
    def __call__(self, img):
        l, a, b = cv2.split(cv2.cvtColor(np.array(img), cv2.COLOR_RGB2LAB))
        return Image.fromarray(cv2.cvtColor(cv2.merge((self.clahe.apply(l), a, b)), cv2.COLOR_LAB2RGB))

class AutoZoom:
    def __init__(self, keep=0.8): self.keep = keep
    def __call__(self, img):
        w, h = img.size
        nw, nh = int(w*self.keep), int(h*self.keep)
        return img.crop(((w-nw)//2, (h-nh)//2, (w+nw)//2, (h+nh)//2))

def get_transforms(sz):
    base = [ApplyCLAHE(3.0), AutoZoom(cfg.ZOOM_FACTOR)]
    return T.Compose(base + [T.RandomResizedCrop(sz, scale=(0.7, 1.0)), T.RandomHorizontalFlip(), T.RandomVerticalFlip(), T.ToTensor(), T.Normalize([0.485, 0.456, 0.
        T.Compose(base + [T.Resize((sz, sz)), T.CenterCrop(sz), T.ToTensor(), T.Normalize([0.485, 0.456, 0.406], [0.229, 0.224, 0.225])])

# 5. Dataset & EarlyStopping Class
class CervicalDataset(Dataset):
    def __init__(self, df, t=None): self.df = df; self.t = t
    def __len__(self): return len(self.df)
    def __getitem__(self, i):
        return self.t(Image.open(self.df.iloc[i]['path']).convert("RGB")) if self.t else Image.open(self.df.iloc[i]['path']).convert("RGB"), torch.tensor(self.df.ilo

class EarlyStopping:
    def __init__(self, patience=5, delta=0):
        self.patience = patience
        self.delta = delta
        self.best_loss = None
        self.counter = 0
        self.early_stop = False

    def __call__(self, val_loss):
        if self.best_loss is None:
            self.best_loss = val_loss
        elif val_loss > self.best_loss - self.delta:
            self.counter += 1
            if self.counter >= self.patience:
                self.early_stop = True
        else:
            self.best_loss = val_loss
            self.counter = 0

# 6. Train/Val Functions
def train_one_epoch(model, loader, criterion, optimizer):
    model.train()
    # Freeze BN to avoid destroying stats with small batches
    for m in model.modules():
        if isinstance(m, nn.BatchNorm2d): m.eval()

    losses = []
    for x, y in loader:
        x, y = x.to(device), y.to(device)
        optimizer.zero_grad()
        loss = criterion(model(x), y)
        loss.backward()
        optimizer.step()
        losses.append(loss.item())
    return np.mean(losses)

@torch.no_grad()
def validate(model, loader, criterion):
    model.eval()
    losses, preds, targets, probs = [], [], [], []

    for x, y in loader:
        x, y = x.to(device), y.to(device)
        logits = model(x)
        loss = criterion(logits, y)

        losses.append(loss.item())
        preds.extend(logits.argmax(1).cpu().numpy())

```

```

        targets.extend(y.cpu().numpy())
        probs.extend(torch.softmax(logits, dim=1).cpu().numpy())

# Calculate Metrics
loss = np.mean(losses)
acc = np.mean(np.array(preds) == np.array(targets))
f1 = f1_score(targets, preds, average='macro')

try: auc = roc_auc_score(targets, np.array(probs), multi_class='ovr', average='weighted')
except: auc = 0.5

return loss, acc, f1, auc, preds, targets

# 7. Main
def main():
    train_df, val_df, test_df = prepare_dataframes()
    train_tfms, val_tfms = get_transforms(cfg.IMAGE_SIZE)

    # Weighted Sampler
    counts = train_df['class_num'].value_counts().sort_index()
    weights = [1.0/counts[l] for l in train_df['class_num']]
    sampler = WeightedRandomSampler(weights, len(weights), replacement=True)

    loaders = {
        'train': DataLoader(CervicalDataset(train_df, train_tfms), batch_size=cfg.BATCH_SIZE, sampler=sampler, num_workers=2),
        'val': DataLoader(CervicalDataset(val_df, val_tfms), batch_size=cfg.BATCH_SIZE, shuffle=False, num_workers=2),
        'test': DataLoader(CervicalDataset(test_df, val_tfms), batch_size=cfg.BATCH_SIZE, shuffle=False, num_workers=2)
    }

    # Model Setup
    model = timm.create_model(cfg.MODEL_NAME, pretrained=True, num_classes=cfg.NUM_CLASSES, drop_path_rate=cfg.DROP_PATH_RATE)
    for p in model.parameters(): p.requires_grad = False

    # Unfreeze Parts
    for p in model.classifier.parameters(): p.requires_grad = True
    for p in model.blocks[-1].parameters(): p.requires_grad = True
    for p in model.blocks[-2].parameters(): p.requires_grad = True

    # Add Dropout to Head
    model.classifier = nn.Sequential(nn.Dropout(p=cfg.DROPOUT_RATE), nn.Linear(model.classifier.in_features, cfg.NUM_CLASSES))
    model = model.to(device)

    criterion = nn.CrossEntropyLoss(label_smoothing=0.1)
    optimizer = torch.optim.AdamW(model.parameters(), lr=cfg.LR, weight_decay=0.01)
    scheduler = torch.optim.lr_scheduler.CosineAnnealingLR(optimizer, T_max=cfg.EPOCHS)

    # Early Stopping Init
    early_stopper = EarlyStopping(patience=cfg.PATIENCE, delta=0.001)

    out_dir = Path(cfg.OUT_DIR) / f"exp_{time.strftime('%m%d_%H%M')}"
    out_dir.mkdir(parents=True, exist_ok=True)

    best_auc = -1.0
    history = {'train_loss': [], 'val_loss': [], 'val_auc': []}

    print(f"\nStarting Training (Max {cfg.EPOCHS} Epochs, Patience {cfg.PATIENCE})...")

    for epoch in range(1, cfg.EPOCHS + 1):
        t_loss = train_one_epoch(model, loaders['train'], criterion, optimizer)
        v_loss, v_acc, v_f1, v_auc, _, _ = validate(model, loaders['val'], criterion) # <--- FIXED: Unpacking 6 items
        scheduler.step()

        history['train_loss'].append(t_loss)
        history['val_loss'].append(v_loss)
        history['val_auc'].append(v_auc)

        print(f"Epoch {epoch:02d} | Loss: {t_loss:.4f} | Val Loss: {v_loss:.4f} | Val Acc: {v_acc:.3f} | Val AUC: {v_auc:.3f}")

        # Save Best
        if v_auc > best_auc:
            best_auc = v_auc
            torch.save(model.state_dict(), out_dir / "best_model.pth")

        # Check Early Stopping
        early_stopper(v_loss)
        if early_stopper.early_stop:
            print(f"Early stopping triggered at Epoch {epoch}!")
            break

    # --- PLOTTING ---
    plt.figure(figsize=(12, 5))
    plt.subplot(1, 2, 1)
    plt.plot(history['train_loss'], label='Train Loss')
    plt.plot(history['val_loss'], label='Val Loss')
    plt.title('Loss')
    plt.legend()

    plt.subplot(1, 2, 2)
    plt.plot(history['val_auc'], label='Validation AUC', color='green')
    plt.title('AUC')
    plt.legend()

    plt.tight_layout()
    plt.savefig(out_dir / "training_plots.png")
    plt.show()

    # --- FINAL TEST ---
    print("\nLoading Best Model for Final Test...")
    model.load_state_dict(torch.load(out_dir / "best_model.pth"))
    test_acc, test_f1, test_auc, preds, targets = validate(model, loaders['test'], criterion) # <--- FIXED: Unpacking

```

```

_, test_acc, test_auc, test_loss, test_f1, test_support = evaluate_model(model, loader, device, patience)

print(f"TEST RESULTS -> Acc: {test_acc:.3f} | AUC: {test_auc:.3f}")
print(classification_report(targets, preds, target_names=["Normal", "Precancer", "Cancer"]))

cm = confusion_matrix(targets, preds)
plt.figure(figsize=(6, 5))
sns.heatmap(cm, annot=True, fmt='d', cmap='Blues', xticklabels=["Normal", "Precancer", "Cancer"], yticklabels=["Normal", "Precancer", "Cancer"])
plt.title('Confusion Matrix')
plt.show()

torch.save(model, "/content/drive/MyDrive/MLProject/final_model_full.pkl")
print("Model Saved.")

if __name__ == "__main__":
    main()


```

Drive already mounted at /content/drive; to attempt to forcibly remount, call drive.mount("/content/drive", force_remount=True).
Device: cuda

Starting Training (Max 25 Epochs, Patience 5)...

Epoch	Loss	Val Loss	Val Acc	Val AUC
Epoch 01	0.8976	0.9295	0.536	0.848
Epoch 02	0.6590	1.0329	0.513	0.854
Epoch 03	0.5697	0.9828	0.557	0.875
Epoch 04	0.5218	1.0100	0.586	0.875
Epoch 05	0.4815	1.0860	0.582	0.869
Epoch 06	0.4526	1.0887	0.604	0.871


Early stopping triggered at Epoch 6!



Loading Best Model for Final Test...

TEST RESULTS -> Acc: 0.636 | AUC: 0.722

	precision	recall	f1-score	support
Normal	0.74	0.76	0.75	34
Precancer	0.38	0.40	0.39	15
Cancer	0.75	0.50	0.60	6
accuracy			0.64	55
macro avg	0.62	0.55	0.58	55
weighted avg	0.64	0.64	0.64	55



Model Saved.

```

# =====
# NOTEBOOK 3: HIGH-CONFIDENCE EVALUATION (Fixed for PyTorch 2.6+)
# =====

```

```

import os
import re
import torch
import numpy as np
import pandas as pd

```

```

from PIL import Image
from pathlib import Path
from torch.utils.data import Dataset, DataLoader
import torchvision.transforms as T
import cv2
import matplotlib.pyplot as plt
import seaborn as sns
from sklearn.metrics import classification_report, confusion_matrix
from google.colab import drive

# 1. Setup
drive.mount('/content/drive')

# --- CONFIGURATION ---
MODEL_PATH = "/content/drive/MyDrive/MLProject/final_model_full.pkl"
TEST_DIR = "/content/via_dataset_ready/test"
CSV_PATH = "/content/drive/MyDrive/MLProject/labels_combine.csv"
CONFIDENCE_THRESHOLD = 0.75 # 75%
IMAGE_SIZE = 384
ZOOM_FACTOR = 0.8

# --- CLASSES ---
class AutoZoom(object):
    def __init__(self, keep_percent=0.65):
        self.keep_percent = keep_percent
    def __call__(self, img):
        w, h = img.size
        new_w = int(w * self.keep_percent)
        new_h = int(h * self.keep_percent)
        left = (w - new_w)/2; top = (h - new_h)/2
        right = (w + new_w)/2; bottom = (h + new_h)/2
        return img.crop((left, top, right, bottom))

class ApplyCLAHE(object):
    def __init__(self, clip_limit=2.0, tile_grid_size=(8, 8)):
        self.clahe = cv2.createCLAHE(clipLimit=clip_limit, tileGridSize=tile_grid_size)
    def __call__(self, img):
        img_np = np.array(img)
        lab = cv2.cvtColor(img_np, cv2.COLOR_RGB2LAB)
        l, a, b = cv2.split(lab)
        cl = self.clahe.apply(l)
        limg = cv2.merge((cl, a, b))
        final = cv2.cvtColor(limg, cv2.COLOR_LAB2RGB)
        return Image.fromarray(final)

class EvalDataset(Dataset):
    def __init__(self, df, transform=None):
        self.df = df
        self.transform = transform
    def __len__(self): return len(self.df)
    def __getitem__(self, idx):
        row = self.df.iloc[idx]
        img = Image.open(row['path']).convert("RGB")
        if self.transform: img = self.transform(img)
        return img, row['class_num'], row['patient_id']

# --- MAIN EVALUATION ---
def evaluate_model():
    device = torch.device("cuda" if torch.cuda.is_available() else "cpu")
    print(f"Using device: {device}")

    # 1. Load Dataframe (Test Set Only)
    master_df = pd.read_csv(CSV_PATH)
    label_map = dict(zip(master_df['filename'], master_df['label']))

    data = []
    if not os.path.exists(TEST_DIR):
        print("Error: Test directory not found. Please unzip data first.")
        # Optional: Uncomment if needed
        !unzip -q "/content/drive/MyDrive/ML Project/Cleaned_Dataset_Processed.zip" -d "/content/via_dataset_ready"
        return

    print("Scanning Test Folder...")
    for f in os.listdir(TEST_DIR):
        if not f.lower().endswith('.jpg'): continue

        # Regex to match augmented naming convention
        clean_name = re.sub(r'_aug_.*', '', Path(f).stem)
        original_key = f"{clean_name}.jpg"

        # Get Patient ID for reporting
        if "IARC" in clean_name:
            match = re.search(r"([0-9]+[A-Z]+|[A-Z]+[0-9]+)", clean_name)
            patient_id = match.group(0) if match else clean_name
        elif "sample" in clean_name:
            patient_id = clean_name.split('_')[0]
        else:
            patient_id = clean_name

        if original_key in label_map:
            data.append({
                "path": os.path.join(TEST_DIR, f),
                "label": label_map[original_key],
                "patient_id": patient_id
            })

    test_df = pd.DataFrame(data)
    mapping = {"Normal": 0, "Precancerous": 1, "Cancerous": 2}
    test_df['class_num'] = test_df['label'].apply(lambda x: mapping.get(str(x).capitalize(), -1))

```

```

test_df = test_df[test_df['class_num'] != -1]

print(f"Found {len(test_df)} test images.")

# 2. Transforms (Must match Training exactly)
eval_tfms = T.Compose([
    ApplyCLAHE(clip_limit=3.0),
    AutoZoom(keep_percent=ZOOM_FACTOR),
    T.Resize((IMAGE_SIZE, IMAGE_SIZE)),
    T.CenterCrop(IMAGE_SIZE),
    T.ToTensor(),
    T.Normalize(mean=(0.485, 0.456, 0.406), std=(0.229, 0.224, 0.225)),
])

test_ds = EvalDataset(test_df, transform=eval_tfms)
test_loader = DataLoader(test_ds, batch_size=16, shuffle=False)

# 3. Load Model with Security Flag Disabled
print(f"Loading model from {MODEL_PATH}...")

# --- FIX IS HERE ---
model = torch.load(MODEL_PATH, map_location=device, weights_only=False)
# -----

model.eval()

# 4. Inference Loop
high_conf_preds = []
high_conf_targets = []
excluded_cases = []

print("\nRunning Inference...")

with torch.no_grad():
    for imgs, targets, patient_ids in test_loader:
        imgs = imgs.to(device)
        logits = model(imgs)
        probs = torch.softmax(logits, dim=1)

        # Get max probability and predicted class
        max_probs, preds = torch.max(probs, dim=1)

        # CPU conversion
        max_probs = max_probs.cpu().numpy()
        preds = preds.cpu().numpy()
        targets = targets.numpy()

        for i in range(len(preds)):
            pid = patient_ids[i]
            confidence = max_probs[i]
            true_label = targets[i]
            pred_label = preds[i]

            # --- CONFIDENCE FILTER ---
            if confidence >= CONFIDENCE_THRESHOLD:
                high_conf_preds.append(pred_label)
                high_conf_targets.append(true_label)
            else:
                excluded_cases.append({
                    'Patient ID': pid,
                    'Confidence': f"{confidence:.2%}",
                    'Predicted': list(mapping.keys())[pred_label],
                    'True': list(mapping.keys())[true_label]
                })

# 5. Results
accepted_count = len(high_conf_preds)
excluded_count = len(excluded_cases)
total_count = accepted_count + excluded_count

print("-" * 40)
print(f"TOTAL TEST IMAGES: {total_count}")
print(f"ACCEPTED (>= {CONFIDENCE_THRESHOLD*100}%): {accepted_count}")
print(f"EXCLUDED (< {CONFIDENCE_THRESHOLD*100}%): {excluded_count}")
print("-" * 40)

if accepted_count > 0:
    target_names = ["Normal", "Precancer", "Cancer"]
    print("\nClassification Report (High Confidence Only):")
    print(classification_report(high_conf_targets, high_conf_preds, target_names=target_names))

    # Confusion Matrix
    cm = confusion_matrix(high_conf_targets, high_conf_preds)
    plt.figure(figsize=(6, 5))
    sns.heatmap(cm, annot=True, fmt='d', cmap='Greens',
                xticklabels=target_names, yticklabels=target_names)
    plt.title(f'Confusion Matrix (Confidence >= {CONFIDENCE_THRESHOLD})')
    plt.xlabel('Predicted')
    plt.ylabel('True')
    plt.show()
else:
    print("No images met the confidence threshold!")

# 6. List Excluded Cases
if excluded_cases:
    print("\n--- EXCLUDED CASES (Low Confidence) ---")
    excluded_df = pd.DataFrame(excluded_cases)
    # Sort by ID for readability
    excluded_df = excluded_df.sort_values('Patient ID')

```

```
print(excluded_df.to_string(index=False))

# Optional: Save to CSV
excluded_df.to_csv("excluded_low_confidence_cases.csv", index=False)
print("\nSaved list to 'excluded_low_confidence_cases.csv'")

if __name__ == "__main__":
    evaluate_model()
```


Drive already mounted at /content/drive; to attempt to forcibly remount, call drive.mount("/content/drive", force_remount=True).
Using device: cuda
Scanning Test Folder...
Found 55 test images.
Loading model from /content/drive/MyDrive/MLProject/final_model_full.pkl...

Running Inference...

TOTAL TEST IMAGES: 55
ACCEPTED (>= 75.0%): 34
EXCLUDED (< 75.0%): 21

Classification Report (High Confidence Only):

	precision	recall	f1-score	support
Normal	0.75	0.96	0.84	25
Precancer	0.00	0.00	0.00	7
Cancer	1.00	0.50	0.67	2
accuracy			0.74	34
macro avg	0.58	0.49	0.50	34
weighted avg	0.61	0.74	0.66	34



--- EXCLUDED CASES (Low Confidence) ---

Patient ID	Confidence	Predicted	True
AAH1	47.91%	Cancerous	Cancerous
AAW1	51.52%	Normal	Normal
ABV1	40.55%	Normal	Precancerous
ACN1	56.37%	Precancerous	Normal
ACQ1	69.70%	Normal	Precancerous
ADB1	74.51%	Normal	Precancerous
ADD1	72.68%	Normal	Precancerous
ADM1	58.89%	Normal	Precancerous
AEH1	63.47%	Normal	Precancerous
AEM1	56.14%	Normal	Cancerous
AFA1	68.29%	Cancerous	Cancerous
AFH1	74.58%	Precancerous	Normal
AGM1	74.32%	Normal	Normal
AHD1	54.59%	Precancerous	Normal
AI01	74.49%	Precancerous	Normal
AJW1	70.25%	Normal	Normal
ANV1	60.90%	Normal	Normal
sample25	70.73%	Precancerous	Cancerous
sample3	57.64%	Normal	Normal
sample45	66.69%	Normal	Precancerous
sample73	54.66%	Normal	Precancerous

Saved list to 'excluded_low_confidence_cases.csv'

```
import os
import re
import math
import pandas as pd
import matplotlib.pyplot as plt
from PIL import Image
from pathlib import Path

# --- CONFIGURATION ---
TEST_DIR = "/content/via_dataset_ready/test"
CSV_PATH = "excluded_low_confidence_cases.csv"
ZOOM_FACTOR = 0.65

class AutoZoom(object):
    def __init__(self, keep_percent=0.65):
        self.keep_percent = keep_percent
    def __call__(self, img):
        w, h = img.size
```

```

        new_w = int(w * self.keep_percent)
        new_h = int(h * self.keep_percent)
        left = (w - new_w)/2; top = (h - new_h)/2
        right = (w + new_w)/2; bottom = (h + new_h)/2
        return img.crop((left, top, right, bottom))

def display_grouped_errors():
    # 1. Load Data
    if not os.path.exists(CSV_PATH):
        print(f"Error: '{CSV_PATH}' not found.")
        return

    df = pd.read_csv(CSV_PATH)
    if len(df) == 0:
        print("No excluded cases to display.")
        return

    # 2. Map Files
    id_to_path = {}
    for f in os.listdir(TEST_DIR):
        if not f.lower().endswith('.jpg'): continue
        clean_name = re.sub(r'_aug.*', '', Path(f).stem)

        # ID Extraction Logic
        if "IARC" in clean_name:
            match = re.search(r"([0-9]+[A-Z]+[A-Z]+[0-9]+)", clean_name)
            pid = match.group(0) if match else clean_name
        elif "sample" in clean_name:
            pid = clean_name.split('_')[0]
        else:
            pid = clean_name

        if pid in df['Patient ID'].values:
            id_to_path[pid] = os.path.join(TEST_DIR, f)

    # 3. Group by (True Label -> Predicted Label)
    groups = {}
    for _, row in df.iterrows():
        key = (row['True'], row['Predicted'])
        if key not in groups:
            groups[key] = []
        groups[key].append(row)

    # 4. Display Logic
    zoomer = AutoZoom(keep_percent=ZOOM_FACTOR)

    # Sort keys for consistent order
    sorted_keys = sorted(groups.keys())

    for true_lbl, pred_lbl in sorted_keys:
        rows_data = groups[(true_lbl, pred_lbl)]

        # Header for the Group
        print("\n" + "="*60)
        print(f"GROUP: True '{true_lbl}' --> Predicted '{pred_lbl}'")
        print(f"Count: {len(rows_data)}")
        print("="*60)

        # Setup Grid
        num_imgs = len(rows_data)
        cols = 4
        rows = math.ceil(num_imgs / cols)

        fig, axes = plt.subplots(rows, cols, figsize=(16, 4 * rows))
        if num_imgs == 1: axes = [axes] # Handle single case
        else: axes = axes.flatten()

        # Plot each image in the group
        for i, row in enumerate(rows_data):
            pid = row['Patient ID']
            conf = row['Confidence']

            if pid in id_to_path:
                try:
                    img = Image.open(id_to_path[pid])
                    img = zoomer(img)

                    # Set color: Green if they match (just low conf), Red if error
                    color = 'green' if true_lbl == pred_lbl else 'red'

                    ax = axes[i]
                    ax.imshow(img)
                    ax.set_title(f"ID: {pid}\nConf: {conf}", fontsize=11, fontweight='bold', color=color)
                    ax.axis('off')
                except Exception as e:
                    print(f"Error loading {pid}: {e}")
            else:
                # Placeholder for missing file
                axes[i].text(0.5, 0.5, "File Not Found", ha='center')
                axes[i].axis('off')

        # Clear remaining empty axes in the grid
        for j in range(i + 1, len(axes)):
            axes[j].axis('off')

        plt.tight_layout()
        plt.show()

display_grouped_errors()

```


GROUP: True 'Cancerous' --> Predicted 'Cancerous'
Count: 2

ID: AAH1
Conf: 47.91%




ID: AFA1
Conf: 68.29%




GROUP: True 'Cancerous' --> Predicted 'Normal'
Count: 1

Error loading AEM1: 'numpy.ndarray' object has no attribute 'imshow'



GROUP: True 'Cancerous' --> Predicted 'Precancerous'
Count: 1

Error loading sample25: 'numpy.ndarray' object has no attribute 'imshow'



GROUP: True 'Normal' --> Predicted 'Normal'
Count: 5

ID: AAW1
Conf: 51.52%



ID: AGM1
Conf: 74.32%



ID: AJW1
Conf: 70.25%



ID: AJP1
Conf: 70.25%



Start coding or generate with AI.



ID: sample3
Conf: 57.64%



